

CIRCULATING & UNSTIRRED WATER BATH

DSB-D & WB-D series



DSB-500D

WB-1000D

WB-500D

Features:

- DSB-500D & DSB-1000D are circulating water bath which have Immersion pump.
- WB-500D & WB-1000D are unstirred water bath.
- Excellent temperature control.
- Large bath opening lets you immerse samples conveniently.
- Inside stainless steel tank.
- Long life immersion stainless pipe heater.
- Quality certificate: ISO 9001, ISO 13485.

Specifications :

Model	DSB-500D	DSB-1000D	WB-500D	WB-1000D
Stirring type	Circulating		Unstirred	
Capacity	10 Liters	20 Liters	10 Liters	20 Liter
Inner tank dimension (mm) (WxDxH)	240x300x150	290x490x150	240x300x150	290x490x150
Temperature setting / readout	Digital / Digital		Digital / Digital	
Temperature control	PID		PID	
Temperature sensor	PT 100Ω		PT 100Ω	
Stability	± 0.2 °C		± 0.2 °C	
Working temperature range	Ambient + 5 ~ 100 °C			
Power watt	1200W		650W	1200W
Optional accessory	Flat lid, gable cover, external circulation equipment, test tube rack		Flat lid, gable cover, test tube rack	
Dimension (mm) (WxDxH)	275x335x190	340x540x190	290x400x240	340x600x240
Weight (kg) (N.W./G.W.)	8.0 / 9.5	10.0 / 11.5	7.0 / 8.5	11.0 / 12.5

Power supply: AC 110V 60Hz (220V 50/60Hz available)



Flat lid



Gable cover



Control head with external circulation equipment

Specifications for test tube racks :

Model	Capacity	Dimension (mm) (WxDxH)
TR-18020	Φ18mm x 20place	103x87x92
TR-16050	Φ16mm x 50place	184x98x71
TR-18050	Φ18mm x 50place	240x107x91
TR-22050	Φ22mm x 50place	244x127x91
TR-18100	Φ18mm x 100place	250x220x92



No.1, Lane 393, Fu De 1 Road, Xizhi Dist., New Taipei City 22149, Taiwan
 T:886-2-26942339. 26945547. 26941738 F:886-2-26941695
 ✉ digilab@ms23.hinet.net www.digisystemlab.com

We reserve the rights to change the specifications without notice.