

dTrite Electronic Digital Titrator

dTrite is a titrator coupled with magnetic stirrer for efficient operation and accurate titration results. It is a high precision, electronically controlled bottle top titrator.

Features

- It works on positive displacement principle.
- Adjustable titration velocity in max. of 16 stages.
- High accuracy of 0.2% with a precision error of 0.07%.
- Continuous display of dispensed volume through remote control panel prevents manual errors.



Specifications

Volume Range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10 μ L
Volume Accuracy	R=0.2% CV=0.07%
Velocity	16 stages
Operating Temperature Range	10°C-30°C
Quality Control	DIN EN ISO 8655
Control Type	External control, including the stirrer and the burette program