

Levo
Pipette Controller

Features

- Designed for simple and precise pipetting control.
- Integrated with 3µm replaceable filter placed within silicone adapter.
- Aspiration and dispensing can be controlled with a sensitive lever.
- Light in weight and avoids the risk of strain.
- Easy for both left and right handed operations.



Compress the suction bellows before attaching the pipette

Button to expel residual liquid when using blow-out pipette

Aspiration and dispensing comfortably achieved with a sensitive lever



3µm replaceable Hydrophobic filter

Silicone adapter for 0.1-100mL pipette tip, fully autoclavable

