

BRAND: PULISITE @ **HUAZHI**

MANUFACTURE: HUAZHI (FUJIAN) ELECTRONIC TECHNOLOGY CO.,LTD

Focusing on the development and manufacture of excellent balances
and laboratory instruments for years




Item No.	Capacity	Readability	Repeatability	Linearity	Pan Size	Housing(mm)	Calibration Way
PTX-FA220S	220g	0.1mg	±0.1mg	±0.2mg	Φ90mm	345×223×331	Internal Cal

PTX Intelligent Internal Calibration Analytical Balance

PTX series using new numeric mainboard: multiple mainboard technology, higher integration project. PTX series balance has many different application, multiple intelligent anti-shock technology, temperature compensation with whole range, multipoint linearity calibration and so on. Highly compatible software program widely satisfy user's advanced request of high precision balance.

■ Main Characteristic:

- ◆ super Monoblock sensor
- ◆ One-key touch internal calibration
- ◆ Six level anti-vibration adjustable function
- ◆ Aluminum down base+ABS top shell
- ◆ Underhook weighing device
- ◆ Weighing results can be input into text file
- ◆ VATN large LCD display
- ◆ Weighing stableness time adjustable
- ◆ Built-in dynamic temperature compensation can fix the impact of weighing result caused by temperature change.
- ◆ 20 units shield: g, mg, ct, oz, ozl, dwt, GN, lb, N, dr, tIT, tIs, tIH, T(Tola), T/A/R(tola/anna/rati), /A/R (tola/Mna/rati), ms (Mesghal), bat, /lb (Parts per pound), kg
- ◆ RS232/USB port can be connected with computer and printer at the same time.
- ◆ Default Weighing mode lock automatically
- ◆ Date and time setting.
- ◆ With display temperature



Super Monoblock Sensor

■ Application Mode

- ✧ Percentage Weighing
- ✧ Dynamic Weighing
- ✧ Counting Function
- ✧ Computing Function
- ✧ High Low Limit Alarm
- ✧ Gross/Net/Tare Weight Weighing
- ✧ Peak Holding Function
- ✧ Density Weighing

