

Touch Screen Constant-Temperature Incubator


Features:

- * Large LCD touch screen.
- * D type: with inner glass door.
G type: with viewing window.
- * Microprocessor controller.
- * Stainless steel chamber.

Safety Functions:

- * Over temperature alarm.
- * Independent limit temperature alarm system.
- * Leakage current and over-voltage protector.

Standard Accessories:



Test port

Optional Accessories:



Base stand



Universal casters



UV lamp



Double-layer Door (D type)



Glass Viewing Window Door (G type)

Technical Parameters:

| Model | BJPX-H54BK(D/G) | BJPX-H88BK(D/G) | BJPX-H160BK(D/G) |
|----------------------------|--|--|---|
| Capacity | 54L | 88L | 160L |
| Temperature Control Method | PT100 | | |
| Temperature Range | RT+5~85°C | | |
| Temperature Stability | ±0.5°C (at 37°C) | | |
| Temperature Uniformity | ±1°C (at 37°C) | | |
| Ambient Temperature | 5~30°C, recommend 25°C±2°C | | |
| Timing Range | 0~999h | | |
| Shelves | 2pcs/adjustable | | |
| Power Consumption | 345W | | 495W |
| Power Supply | AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional) | | |
| Internal Size(W*D*H)mm | 395*346*395 | 500*364*500 | 558*518*600 |
| External Size(W*D*H)mm | D type: 536*588*681 G type: 522*626*587 | D type: 650*573*793 G type: 627*621*692 | D type: 718*732*890 G type: 685*779*790 |
| Package Size(W*D*H)mm | D type: 606*648*818 G type: 592*676*724 | D type: 720*630*930 G type: 697*671*829 | D type: 788*792*1027 G type: 755*829*927 |
| Gross Weight(kg) | D type: 51 G type: 44 | D type: 65 G type: 51 | D type: 89 G type: 69 |

Biochemistry Incubator

Features:


- * LCD touch screen.
- * Double door design.
- * Polished stainless steel chamber.
- * Round angle structure, easy to clean.
- * Universal castors for easy movement.
- * CFC-free environmental friendly refrigerant.
- * Auto defrost and manual defrost both available
- * Microprocessor PID temperature control: accurate and reliable.
- * BOD socket integrated inside for additional small equipment.

Safety Functions:

- * Temperature protection threshold value can be set.
- * Auto-off protection and alarm for abnormal situation.
- * Parameters can be saved automatically when power off or system error.



LCD touch screen



Built-in printer



Test port



Double door

Technical Parameters:

| Model | BJPX-B100 | BJPX-B150 | BJPX-B200 | BJPX-B250 |
|----------------------|--|--------------|--------------|--------------|
| Capacity | 100L | 150L | 200L | 250L |
| Temp. Range | 0~60°C | | | |
| Temp. Fluctuation | ±1°C | | | |
| Temp. Uniformity | ±2°C at 37°C | | | |
| Shelves | 2 pcs/ adjustable | | | |
| Time Range | 0~99h | | | |
| Refrigerant | R134a (CFC free, eco-friendly) | | | |
| Working Type | Continuous | | | |
| Standard Accessory | BOD socket, built-in printer, 2 test ports, lock | | | |
| Optional accessories | USB port | | | |
| Consumption | 1200W | | 1400W | 1500W |
| Power Supply | AC220V±10%, 50/60Hz(Standard); AC110±10%, 60Hz(Optional) | | | |
| Internal Size(mm) | 500*350*580 | 550*401*690 | 600*450*750 | 660*450*850 |
| External Size(mm) | 644*647*1170 | 690*692*1270 | 742*743*1340 | 802*743*1440 |
| Packing Size(mm) | 769*772*1327 | 815*817*1327 | 865*867*1497 | 865*927*1597 |
| Gross Weight(kg) | 105 | 120 | 130 | 145 |