


JY-SPCT Horizontal Electrophoresis Tank



- Abundant buffer not only guarantees cooling effects, but also keeps the pH value stable during the whole experiment process.
- Removable electrodes enable maintenance and replacement easy and quick.
- With background color design on the tray, convenient to adding samples.
- It's more convenient to open the lid with the easy-open buttons.
- It prevents buffer volatilizing and electricity leakage with transparent lid.
- Auto-switch-off when the lid is moved.

- With the drift prevention function on the tray, effectively prevent the band bending.

JY-SPCT Technical Specifications

Gel size (W×L)	120×120 (mm)
	120×60 (mm)
	60×120 (mm)
	60×60 (mm)
No. of samples	1~100
Total buffer volume	550 (ml)
Dimension (L×W×H)	300×160×75 (mm)
Weight	1.1 (kg)

JY-SPCT Packing List

Parts name	Quantity
Lid with cables	1
Body tank with electrodes	1
120×120mm Gel tray	1
120×60mm Gel tray	1
60×120mm Gel tray	1
60×60mm Gel tray	2
Gel box	1
1.5mm thickness 1/2/3 wells comb	1
1.0mm thickness 6/13 wells comb	2
1.0mm thickness 8/18 wells comb	2
1.0mm thickness 11/25 wells comb	2
level adjusting button	4