




Ultrasonic Cleaner DR-MH Series

Derui International Limited  Ultrasonic Cleaner	Contact & info.	Tel: 06-6114-0543 Email: provsci9@gmail.com
Customer		

Approval /Certificate: CE & RoHS & FCC

Model	Photo	ACCESSORIES
DR-MH Series		

Feature	1)The stainless steel tank used and the reliable technology prolong the service life of the units 2)With mechanical timer & heater control 3) Ultrasonic Timer 0~60 Minute, infinity 4)Temperature adjust from 30°C to 80°C 5) High-performance transducer system
----------------	---

Specifications

Model No.	Ultrasonic power	Heating power	Ultrasonic frequency	Drain	Tank Vol/Size(mm) LxWxD	TransducerPC	Unit Weight	Packing	Packing size(cm)
DR-MH07	60 W	60 W	40 KHz	NO	0.7/150x86x65	1 x 60w	2.00KG	12 unit/ctn	57x50x54
DR-MH13	60 W	60 W	40 KHz	NO	1.3L/150x138x65	1 x 60w	2.14KG	4 unit/ctn	53x53x35
DR-MH20	60 W	100 W	40 KHz	NO	2L/150x138x100	1 x 60w	2.18KG	4 unit/ctn	53x53x35
DR-MH22	120W	100 W	40 KHz	NO	2.2L/238x138x65	2 x 60w	2.50KG	4 unit/ctn	49x34x60
DR-MH30	120 W	150W	40 KHz	NO	3L/238x138x100	2 x 60w	3.20KG	4 unit/ctn	49x34x60
DR-MH40	120 W	220W	40 KHz	NO	4L/302x152x100	2 x 60w	3.75KG	2 unit/ctn	43x57x34
DR-MH45	120 W	220W	40 KHz	NO	4.5L/238x138x150	2 x 60w	4.50KG	2 unit/ctn	49x34x39
DR-MH60	180 W	270W	40 KHz	NO	6L/302x152x150	3 x 60w	4.72KG	2 unit/ctn	57x42x39
DR-MH100	240 W	310W	40 KHz	YES	10L/302x239x150	4 x 60w	5.20KG	1 unit/ctn	47x37x43
DR-MH120	300 W	380w	40 KHz	YES	12L/490x138x150	5 x 60w	5.80KG	1 unit/ctn	65x25x47
DR-MH130	300 W	420W	40 KHz	YES	13L/302x239x200	5 x 60w	6.50KG	1 unit/ctn	46x36x47
DR-MH150	360 W	470W	40 KHz	YES	15L/325x295x150	6 x 60w	7.50KG	1 unit/ctn	47x42x47
DR-MH220	360 W	470W	40 KHz	YES	22L/502x298x150	6 x 60w	8.75KG	1 unit/ctn	65x42x39
DR-MH280	480 W	620W	40 KHz	YES	28L/502x298x200	8 x 60w	10.00KG	1 unit/ctn	65x42x39