

GENERAL LAB

PROMOTION

GENERAL LAB

ลดแล้ว
ลดอีก



Vacuum Pump

LCH-GMP-15

- Display Mode: Pressure gauge
- Ultimate Pressure (mbar): -0.085
- Positive Pressure (mPa): 0.2
- Pumping Speed (L/min): 15
- Pump Heads: 2

11,500.-



Digital Vortex Mixer

MX-S+

- Speed Rang: 100-3,000 rpm
- Time: 1S-99min59s
- display: LCD

8,000.-



High Speed Mini Centrifuge

D2012S

- Speed Rang: 300-15,000 rpm
- Time: 30S-99min
- Capacity: 1.5/2 ml x 12

Included Rotor 19400158

26,000.-



High Speed Refrigerated Centrifuge

DC1616R

- Speed Rang: 16,000 rpm
- Time: Timing/Hold
- Temperature Setting: -10-40°C
- Capacity: 4 x 400 ml

Included Rotor 19400155

350,000.-



Magnetic Stirrer With Electric Heating Mantle

MS-H280-B

- Speed Rang: 200-1,500 rpm
- Temp Max: 280 °C
- Size: 100,250,500,1000 mL

14,000.-



Balance

PTY-223

- Capacity: 220g
- Readability: 1mg
- Repeatability: ±1mg
- Internal: Calibration

41,000.-



Pipetting,Mixing,Stepper and Dilution

Multi functional Electronic Pipette

dPette+

- 0.5 - 10
- 5-50
- 30-300
- 100-1000

13,000.-



Pipetting,Mixing

Simple Electronic Pipette

dPette

- 0.5 - 10
- 5-50
- 30-300
- 100-1000

10,000.-



science.shop

บริษัท โปรวิชั่น ไซแอนติฟิค จำกัด

252 ถนนสวนสยาม แขวงคันนายาว เขตคันนายาว กรุงเทพฯ 10230
โทร : 0 2517 5772-4 / โทรสาร : 0 2517 5772-4 ต่อ 818
252 Suan Siam Road , Khannayao , Khannayao Bangkok 10230
Tel.0 2517 5772-4 / Fax.0 2517 5772-4 - 818 www.provsci.com



Tel. 0 2517 5772-4



Provision Scientific



provsci9@gmail.com



www.provsci.com

CONTACT

1 July - 31 December 2025

All items: Warranty 1 year/Excluded vat 7%



ตลาดวิทยาศาสตร์

@provsci



ADD FRIEND

1



ONILAB

Orbital Shaker**SK-0180-S**

- Speed: 40-200 rpm
- Max. load: 2kg
- Display: LED

Included platform 1

6,300.-

ONILAB

Benchtop-Cond.Meter**EC100**

- Cond.Range: 0-200.0mS
- TDS: 0-100.0g/L
- Salinity: 0-100ppt
- Resolution: 0.001
- accuracy: $\pm 1.0\%FS \pm \text{one digit}$

20,700.-

LANPHAN

Autoclave**LPLS-50LD**

- Chamber: 50L, $\varnothing 340 \times 550 \text{mm}$
- Temp: 134°C
- Time: 0-99 min

75,600.-

LANPHAN

Autoclave**LPYX-18HDD**

- Chamber volume: 18L, $\varnothing 280 \times 260 \text{mm}$
- Working temp.: 126°C
- Max.safety pressure: 0.165MPa
- Time: 0-99min

27,000.-

BestScope

Microscope**BS-08030B**

- Viewing Head : Inclined Binocular head, 360° Rotatable, Interpupillary distance 55-75mm
- Eyepiece : WF10x/22mm
- Objective : 1x-4x
- Incident Illumination : Halogen lamp 12V/10W Fluorescent 10W
- Transmitted Illumination : Halogen lamp 12V/10W
- Accessories : High quality Jewel clamp : Darkfield working stage (with Iris diagram)

30,500.-

LACHOI

Oven**LAS-VD-90**

- Display Mode: LCD
- Vacuum (pa): 133
- Temperature Control Range (°C): RT+10-280
- Timing Range (min): 0-9999
- Shelves: 3

115,500.-

GONDO EZDO

PL-700 Series Bench Top Meter**PL-700PVS**

- Range: -2.00~16.00 pH
- Accuracy: $\pm 0.01 + 1 \text{ digit}$
- Resolution: 0.01 pH
- Compensation: ATC: 0~100 °C

23,500.-

XingYun

Moisture**XY-100MW-A**

- Capacity : 110g
- Operating Temp. : 5°C-35°C
- Display : High Definition Backlit LCD
- Pan size : 90 mm
- Temp. Sensor : PT100
- Temp. Range : 40°C - 199°C
- Moisture Readability : 0.01%
- Dry Readability : 0.01%
- Halogen lamps / HBM Sensor

42,500.-

LACHOI

Viscometer**LCH-NDJ-8S**

- Display LCD
- Measurement range 1-2x10⁶ mPa.s
- Spindle specification No.1-4 Spindle
- Spindle Speed 0.3/0.6/1.5/3/6/12/30/60 rpm

32,500.-

Lab control

Orbital Shaker

- Temp.: 0 to 50 °C
- Speed: 40-280 rpm
- Volume(L): 100
- Std.Holder: Spring Tray
- Shaking mode: Orbital

115,000.-

BIOBASE

Incubator**BJPX-L400II**

- Capacity : 402L
- Temperature Range: 0~65°C
- Timing Range: 0-9999 min/h
- Temperature control Range : With illumination : 10~50 °C; Without illumination : 5~50 °C
- Timing Range : 0-99 h
- LCD display
- 3-side lighting
- Intensity of Illumination : 0-12000 Lux
- Illumination Adjustment Grade : BJPX-LII series : 5-grade adjustable
- Standard : Shelf*3

230,000.-

BIOBASE

Incubator**BJPX-I-150**

- Capacity: 150L
- Temperature Range: 0~65°C
- Timing Range: 0-9999 min/h
- Shelves : 2pcs/adjustable

57,600.-