

HALLOWEEN CLEARANCE



รับฮาโลวีนกับ โปรโมชั่นปล่อยผี

วันนี้ - 31 ธันวาคม 2025
(หรือจนกว่าสินค้าจะหมด)

Ductless Fume Hood



SCITEK

MODEL **FMH-800DL**
Overall Dimension (W×D×H)
800×620×1800 mm
Working Space (W×D×H)
780×540×775 mm
Maximum Opening Height
650 mm
Airflow Velocity
0.4~0.6 m/s, adjustable

ราคา
69,900.-

Benchtop Centrifuge



HWLAB

MODEL **HW6**
Speed range
1,000 ~12,000 rpm
Timer
15s ~ 99min/continuous operation
Max. RCF
approx: 6,900xg

ราคา
8,500.-



Filtered Storage Cabinet

SAFOD

MODEL **SF-Z030**
Overall Dimension (W×D×H)
600×460×1650 mm

ราคา
89,900.-

Galvanized Shelf:	Fan	Filter
1 pc	1 pc	2 pcs A/C filters

“ตัวเครื่องมาพร้อม”
แผ่นกรองถ่านกัมมันต์



Wavelength
control knob

V-1100 Single-beam Visible Spectrophotometer

Spectral Bandwidth
4.0 nm

Light Source and Detector
H lamp, Si-photodiode

Wavelength Accuracy/Repeatability
±2.0 nm / 1.0 nm

Photometric Mode
T, A, E, C

Overall Dimension (W×D×H)
360×280×190 mm

Wavelength Range
325~1000 nm

Photometric Display Range
-0.3~3.0 A

Photometric Accuracy/Repeatability
±0.5% T / 0.2% T

Stray Light:
≤0.1% T (360 nm NaNO₂)

Macy

ราคา
27,900.-





- ✓ Degas
- ✓ Sweep
- ✗ Pulse
- ✗ Power



- ✗ Degas
- ✗ Sweep
- ✗ Pulse
- ✗ Power

Ultrasonic Cleaner

Model	Dimension (LXWxD)	Tank Capacity	ราคา
DR-DS30	238×138×100 mm	3L	9,000.-
DR-DS40	302×152×100 mm	4L	10,500.-
DR-DS60	302×152×150 mm	6L	15,500.-
DR-DS100	490×138×150 mm	10L	19,500.-
DR-MH30	238×138×100 mm	3L	5,900.-
DR-MH40	302×152×100 mm	4L	6,900.-
DR-MH60	302×152×150 mm	6L	8,900.-
DR-MH120	490×138×150 mm	12L	12,900.-
DR-MH150	325×295×150 mm	13L	14,900.-
DR-MH280	502×298×200 mm	28L	19,900.-



Derui®

Auto Dehumidification Dry Cabinet



AD-80S



AD-100S



AD-300S




AD-600S

Model	Dimension (W×D×H)	Capacity	Shelf	Dehumidifier	ราคา
AD-80S	380×390×616 mm	80 L	3 pcs, plastic	1 pc, <8 W	9,500.-
AD-100S	380×390×770 mm	100 L	3 pcs, plastic	1 pc, <8 W	10,500.-
AD-300S	850×430×965 mm	300 L	3 pcs, steel	2 pcs, <15 W	22,500.-
AD-600S	850×430×1835 mm	600 L	4 pcs, steel	4 pcs, <30 W	45,500.-

ANDBON

Rotary Evaporator EV400Vac

SET



Color Touch Control Unit

- ▶ Bidirectional rotation, 20-300 rpm
- ▶ Motorized lifting, 16 cm stroke
- ▶ Built-in vacuum controller

Vertical-Type Condensation Unit

- ▶ Cooling surface: 1,500 sq cm
- ▶ Manual tilt adjustment: 0° to 60°

Main Unit

- ▶ Epoxy-coated aluminum-steel casting
- ▶ PTFE & FKM seal materials
- ▶ Dimension (mm): L600×W400×H900

Heating Bath (5 L)

- ▶ PID control, ambient-230°C ±1°C
- ▶ Heating power: 1,300 W
- ▶ Heating surface: PTFE-coated SS 304L
- ▶ Outer body: PBT

ชุด EV400VAC ฟรี Water Chiller และ Vacuum Pump



Oil-Free Vacuum Pump

VP18R Plus

- ▶ Pump capacity: 18 L/min @20mbar
- ▶ Motor power: 140 W
- ▶ Dimension (mm): L294×W170×H240

39,900 บาท **FREE**



Recirculating Water Chiller

H50-500

- ▶ PID control, -5°C~35°C ±0.3°C
- ▶ Pump capacity: 3 L/min @10psi
- ▶ Reservoir volume: 1.8 L
- ▶ Dimension (mm): L480×W250×H500

92,500 บาท **FREE**

= ราคาคู่มือ =
315,000.-

LabTech